



Your Practice. Your Preference.
Your Productivity Solution.

High Performance DR

The versatile DR system offers images within seconds and fits easily into smaller x-ray rooms.

Vizion Portable DR

- Weight: 8.6lbs
- 14"x17" Imaging Area
- X-Ray Converter - CsI:TI
- Sensor Type - a-Si TFT - PIN Diode
- X-Ray Voltage Range - 40-150kVp
- A/D Conversion - 14 bits
- Image Data Interface - Gigabit Ethernet
- Power Supply - AC100-120V/220-240V, 50/60Hz

Panel Features

- Unsurpassed Flat Panel Technology
- Superior Image Quality
- High Sensitivity
- Sleek Slim Design

Panel Specs

- Pixel Size - 150x150 μ m
- Active Array - 2400 x 2880
- Active Area - 360x432mm

Software Features

- Window level
- Rotate
- APR
- Repeat Reject
- Archive
- Technique Chart
- Patient ID
- Body Region
- Image Preview
- DICOM Print
- DICOM Store

Computer Specs

- Intel Core 2 Duo E7500/2.93GHz
- Dual Core
- 3 GB RAM 800MHz
- 320GB Standard Hard Drive

Monitor Specs

- Single 1.3MP 19" LCD Monitor
- Resolution 1280 X 1024
- 8000:1 Ratio



Viztek • 2217 US Highway 70 East, Garner, NC 27529
Tel: 800-366-5343 • Email: Info@viztek.net