



Your Practice. Your Preference.
Your Productivity Solution.

VERSATILE SOLUTION FOR SUPERIOR IMAGES

Viztek's Wireless DR is portable and lightweight, ideal for all radiology specialties. The durable detector features a convenient handle for transfer, while strong enough to handle accidents that come with constant movement. With image display within seconds, the DR system offers easy integration for workflow enhancement.

Features and Benefits

- Portable and wireless
- Low power WiFi connection means no interruption with machines nearby
- Also available by cable connection where WiFi is not available
- Allows freedom to position special and difficult angles, but also allows for the flexibility to move from room to room in different departments
- Charges on a wall-mounted docking station
- Cesium iodide scintillator, reduces dose exposure and produces high-quality images
- Less than three seconds from exposure to image availability/presentation – enhance workflow
- Easy integration

Technical Specifications

- Size: 14" x 17"
- Pixel Size: 144 μm
- Pixel Matrix: 2372 x 3000 pixels
- AD Conversion: 16 bits
- Generator Volt Range: 40-150 kVp
- Preview Display: 3 seconds
- High Resolution Image Display: 10 seconds
- Battery Duration: 2 hours
- Wireless Connection: WiFi
- Communication Interface: Ethernet
- DC Input Voltage: 24 V
- Weight: 10.58 lbs
- Connection Cable Length (between detector and station): 22.9 ft

Accessories

- Grid Frame
- Battery Charger – Used with bucky tables or wall-stands for continuous detector charging



6491 Powers Avenue • Jacksonville, FL 32217 • Phone: 800.366.5343 • Fax: 904.448.9936
Email: info@viztek.net • www.viztek.net